



ACRO™ FOB Rapid Test Cassette

A rapid, one step test for the qualitative detection of Human Occult Blood in feces.

- Many diseases can cause hidden blood in the feces. This is also known as **Fecal Occult Blood (FOB)**, **Human Occult Blood**, or Human Hemoglobin.
- In the early stages, gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis, and fissures may not show any visible symptoms, **only occult blood**.



Features & Benefits

- Results in 5 minutes
- Easy to use
- No effect by diet.
- Selectively detect FOB at **50ng/ml** (6µg/g) and **100ng/ml** (10µg/g) or higher in feces.

ACRO™ FOB Rapid Test Cassette (Feces)	Sensitivity	Specificity	Accuracy	Format	Package
FOB 50	95.0%	99.5%	98.6%	Cassette	25T
FOB 100	95.5%	99.7%	98.0%	Cassette	25T