


Marijuana and Morphine

One Step Drug Screen Test (Strip or Device)




Test Procedure


STRIP TEST



Step 1. Immerse the strip into urine

Step 2. Read result at 5 minutes

DEVICE TEST



Step 1. Add 3 full drops of urine

Step 2. Read result at 5 minutes



Ordering Information (CONTINUED)

CE Marked

Product Description	Specimen	Format	Catalog No.	Cut off	Kit Size
Morphine (MOP)	Urine	Strip	DMO-101	300 ng/mL	50
		Device	DMO-102		40
Marijuana (THC)	Urine	Strip	DTH-101	50 ng/mL	50
		Device	DTH-102		40