



ACRO™ Influenza A&B Rapid Test Cassette

is a rapid chromatographic immunoassay for the qualitative detection of Influenza A and B antigens in nasopharyngeal swab, throat swab or nasal aspirate specimens


Features & Benefits

- Results in 15 minutes
- Easy to use
- Sample types:
 - ▶ Nasopharyngeal swab
 - ▶ Throat swab
 - ▶ Nasal aspirate specimen

Benefits of Rapid Diagnosis

- Guide appropriate therapy
- Avoid unnecessary antibiotic use
- Reduce cost of treatment


Direction for use



Read the result at 15 minutes

Note: Do not read result after 20 minute

Interpretation of Results



ACRO™ Influenza A Performance vs. RT-PCR

Sample	Sensitivity	Specificity	Accuracy
Nasopharyngeal swab	99.0%	98.9%	98.9%
Throat swab	95.1%	99.3%	98.1%
Nasal aspirate	100%	99.2%	99.3%

ACRO™ Influenza B Performance vs. RT-PCR

Sample	Sensitivity	Specificity	Accuracy
Nasopharyngeal swab	97.7%	99.0%	98.6%
Throat swab	94.2%	99.4%	97.8%
Nasal aspirate	97.9%	99.4%	98.8%

See more about ACRO™ Influenza A&B performance in Package insert.