



ACROTM Leptospira lgG/lgM Rapid Test Cassette

A rapid test for the qualitative detection of IgG and IgM antibodies to Leptospira interrogans

- Leptospirosis is caused by Leptospira interrogans which spread by rodents, dogs, livestock and wildlife.
- Cases of leptospirosis can increase after hurricanes or floods when people may be exposed to contaminated water or use it for drinking or bathing.



Features & Benefits

Results in 15 minutes

i+MED LABORATORIES CO., LTD.

- Easy to use
- Sample types : fingerstick, whole blood, serum or plasma

ACRO™ Leptospira IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma)	Sensitivity	Specificity	Accuracy	Format	Package
Leptospira IgM	90.0%	98.6%	98.3%	Cassette	25T
Leptospira IgG	93.8%	98.7%	98.3%	Cassette	25T





