


THC (Marijuana) Rapid test (Urine)

Test Procedure


Dipstick



Step 1. Immerse the strip into urine

Step 2. Read result at 5 minutes

Cassette



Step 1. Add 3 full drops of urine

Step 2. Read result at 5 minutes



Ordering Information (CONTINUED)

CE Marked

Product Description	Specimen	Format	Catalog No.	Cut off	Kit Size
Marijuana (THC)	Urine	Dipstick	ACR-DTH-101	50 ng/mL	50
		Cassette	ACR-DTH-102		25