



i+LAB

Methamphetamine Test

Cut-off 1000 ng/ml

RAPID DIAGNOSTIC PRODUCT


DRUGS OF ABUSE

i+LAB

- Take the easy step for methamphetamine detection
- Rapid screening in urine sample
- Reliable test result


TEST PROCEDURE

METHOD FOR DEVICE



METHOD FOR STRIP

Immerse the test strip vertically for 30 seconds Avoiding the level of urine specimen over the maximum line on the test strip



ชุดทดสอบ	เลขที่ใบรับแจ้ง	ขนาดบรรจุ
i+LAB Methamphetamine Test Strip, cut off 1000 ng/ml (ชนิดแผ่นทดสอบ)	จพ. 3/2561	50 Tests/Kit
i+LAB Methamphetamine Test Device, cut off 1000 ng/ml (ชนิดถ้วยทดสอบ)	จพ. 2/2561	25 Tests/Kit 40 Tests/Kit 50 Tests/Kit