

i+LAB Methamphetamine Test Cut-off 500 ng/ml


RAPID DIAGNOSTIC PRODUCT DRUGS OF ABUSE



- Take the easy step for methamphetamine detection
- Rapid screening in urine sample
- Reliable test result


TEST PROCEDURE

METHOD FOR DEVICE



METHOD FOR STRIP

Immerse the test strip vertically for 30 seconds. Avoiding the level of urine specimen over the maximum line on the test strip.



ชุดทดสอบ	เลขที่ใบรับแจ้ง	ขนาดบรรจุ
i+LAB Methamphetamine Test Strip, cut off 500 ng/ml (ชนิดแผ่นทดสอบ)	จพ. 7/2561	50 Tests/Kit
i+LAB Methamphetamine Test Device, cut off 500 ng/ml (ชนิดถ้วยทดสอบ)	จพ. 3/2562	25 Tests/Kit 40 Tests/Kit 50 Tests/Kit