



dFlow

Electronic Bottle Top Dispenser

Features

- Easy to operate, light in weight and compact design.
- Operation through separated control panel minimizes the repetitive strain injury.
- Two modes of operation: Dispenser or Stepper.
- Operational volume ranges from 0.1-99.9mL.

Specifications

Volume range	0.1-99.9mL, max piston lift of 10mL, resolution at 100 μ L
Volume accuracy	Dispenser: R = 10mL \pm 30 μ L; CV = 10ml \pm 10 μ L Stepper: R = 1mL \pm 6 μ L; CV = 1mL \pm 9 μ L
Operating temperature Range	10 $^{\circ}$ C-30 $^{\circ}$ C
Quality standard	DIN EN ISO 8655
Contro type	Separated controller