


Trust
Easy Check

Implant Stability Measuring Device

Any Angle

Disposable Vinyl Cap

Similar Osstell Number



Dentium
For Dentists By Dentists

Vinyl Cap

Prevent cross-infection by using disposable plastic cap (sold separately)



Dentium ISV Scale

	60	65	70	
	Low	Medium	High	→
Indication	Implant Risk Monitor ISV	Full Splint IL	Partial	Single
Surgical Protocol		2-Stage	1 or 2-Stage	1-Stage
Restorative Protocol		Traditional	Early	Immediate



Product Name		Easy Check
Body	Weight(g)	94
	Size(mm)	26 x 26 x 224
Holder	Weight(g)	260
	Size(mm)	79 x 111 x 99
Tolerance		±3
Angled		0~90°
Time		Less than 3s
Heating Strength		Less than 1.0N
Timer		Auto power off if not used for more than 1 minute
Rated volatge		3.7V / 0.5A
Charger		USB C-Type
Battery Capacity		420mA
Storage Condition		Temperature 0~40°C / Humidity 30~85%

EYCK- 2020.12[Rev. 0]