

Serological pipettes are extremely efficient when it comes to mixing cell suspensions. But, you can also use them to mix reagents and chemicals. These pipettes are useful in moving colonies of cells for additional analysis during treatment or isolation of cell cultures. Furthermore, you can use these pipettes to layer reagents to create density gradients like a ficoll gradient that refines blood plasma.

Features:

- Easily identified with color code
- Optical-grade polystyrene (PS)
- Cotton plugged for safety
- 2 graduation mark
- Individual wrapped
- Accury: ±2%
- Sterile, non-Pyrogenic, Dnase/Rnase free

Sterile in Individual wrapped

Order Information:

Size	Packaging	Min. graduation	Color	Cat No.
1mL	300 pcs/box	1/100mL	Yellow	SGI080
2mL	300 pcs/box	1/100mL	Green	SGI081
5mL	200 pcs/box	1/10mL	Blue	SGI082
10mL	200 pcs/box	1/10mL	Orange	SGI083
25mL	100 pcs/box	2/10mL	Red	SGI084
50mL	75 pcs/box	5/10mL	Purple	SGI085

