

No.	Rotor	Max. Capacity	Max. Speed(rpm)	Max. RCF(×g)	Adapter
No.1	Tube Rack	4*16*10ml	4000	2810	/
		4*20*7ml PP tube		3139	/
		4*30*5ml blood collection tube(80mm)		3150	/
		4*30*5ml blood collection tube(106mm)		3150	120*2ml
		4*24*5ml blood collection tube(106mm)		3032	96*2ml
		4*24*5ml blood collection tube(100mm)		3032	96*2ml
		4*20*5ml blood collection tube(106mm)		3050	80*2ml
		4*24*2ml blood collection tube		2585	/
		4*20*2ml blood collection tube(80mm)		3050	/
		4*25*1.5ml		2808	/
No.2	Swing Rotor Support	4*750ml	/	/	/
	Hanging Cup	4*750ml	/	/	/
	Adapter	4*3*100ml	4000	3380	/
		4*2*100ml		3360	/
No.3	Swing Rotor Support	4*250ml	/	/	/
	Hanging Cup	4*250ml	/	/	/
	Adapter	4*1*100ml	4000	2379	/
		4*1*50ml		2379	/
		4*9*10ml		2379	/
	Tube Rack	4*4*50ml	4000	3500	/
		4*2*50ml(▲)		3023	/
		4*20*5ml		3139	/
		4*18*5ml blood collection tube(L106mm)		2808	72*2ml
		4*12*5ml blood collection tube(L75mm)		2808	/
		4*12*5ml blood collection tube(L100mm)		2808	48*2ml
		4*12*5ml blood collection tube(L106mm)		2808	48*2ml
		4*10*5ml blood collection tube(L106mm)		2808	40*2ml
No.4	Swing Rotor Support	SUS304 Support	/	/	/
	Hanging Basket	4*2*100ml	4000	3175	/
		4*1*100ml		3007	/

No.	Rotor	Max. Capacity	Max. Speed(rpm)	Max. RCF(×g)	Adapter
No.4	Hanging Basket	4*2*50ml	4000	3022	50ml(▲),15ml(●)
		4*1*50ml		2846	
		4*8*15ml		3030	15ml(▲),10ml(▲), 10ml blood
		4*6*15ml		2705	collection tube, 5ml blood
		4*4*15ml		2705	collection tube, 2ml blood collection tube(L75)
		4*8*10ml		3030	/
		4*6*10ml		2705	/
		4*4*10ml		2705	/
No.5	Microplate Rotor	2*2*96well	4000	2300	/
No.6	Auto Uncap Swing Rotor	4*12*5ml	4000	3580	/
		4*12*2ml		3580	/
		4*20*5ml		3580	/
		4*20*2ml		3580	/
		4*24*5ml		3580	/
		4*24*2ml		3580	/
No.7	Angle Rotor	12*10ml	5000	3440	12*5ml
No.8	Angle Rotor	18*10ml	5000	3516	18*5ml
No.9	Angle Rotor	24*10ml	5000	2200	24*5ml
No.10	Angle Rotor	8*15ml	5000	3300	/
No.11	Angle Rotor	12*15ml	5000	2200	/
No.12	Angle Rotor	12*20ml	5000	3360	/
No.13	Angle Rotor	4*50ml	5000	3400	/
No.14	Angle Rotor	6*50ml	5000	3400	/
No.15	Angle Rotor	4*100ml	4000	2800	/

Note: The default is round-bottom tube, with mark (▲) is sharp-bottom tube. The swing rotor should select support+hanging cup+adapter or support+tube rack.